


Q&A

What is New about
an innovative dentin desensitizer,
NANOSEAL



nanoseal®

INDEX



Q&A about **Product Characteristics**

Q1. What type of dentin desensitizer is NANOSEAL?	3
Q2. What are the new characteristics of NANOSEAL?	4
Q3. What is the composition and pH of NANOSEAL ? Does it smell?	4
Q4. How is the acid-resistant nanoparticle layer formed?	5
Q5. What is the composition of the acid-resistant nanoparticle layer?	6
Q6. How can resistance to acids of the nanoparticle layer be provided?	6
Q7. How strong is the acid resistance of the nanoparticle layer?	6
Q8. Does NANOSEAL contain fluorine?	6

Q&A about **Usage**

Before application

Q1. How long can the mixture of Liquid A and Liquid B be used?	7
Q2. Is there any problem if the collection ratio of Liquid A and Liquid B is not equal?	7
Q3. Is any pre-application treatment required for the tooth surface other than cleaning?	7

During application

Q4. How should NANOSEAL be applied?	7
Q5. Should the surface be recoated?	8
Q6. Is there any problem if NANOSEAL contacts to the gingiva during application?	8
Q7. Do patients feel a bitter taste during application?	8
Q8. Does the tooth hurt during application?	8

After application

Q9. Can the surface be rinsed with water immediately after the application?	8
---	---

Q&A about **Desensitizing effect**

- Q1. What desensitizing effects of NANOSEAL have been reported against hypersensitivity? 9
- Q2. Will additional application be beneficial if no benefit has been observed with a single application? ... 9
- Q3. I am afraid that the nanoparticle layer may be removed by brushing..... 10
- Q4. Can NANOSEAL be applied to the root surface with the remaining cementum?..... 10
- Q5. Does NANOSEAL have the same effects if used on the surface previously coated with
a resin-based dentin desensitizer?..... 10

Q&A about **Application to the Prepared Tooth**

- Q1. Can NANOSEAL be used for the prepared tooth? 11
- Q2. Will the bond strength of resin cement and bonding materials be decreased on the prepared
tooth surface coated with NANOSEAL? 11
- Q3. Is phosphoric acid etching on the NANOSEAL-coated surface (enamel) effective? 11
- Q4. Will the fitting of dental prostheses be affected by NANOSEAL? 12
- Q5. Can temporary seal and cementation materials be used after the application of NANOSEAL? 12
- Q6. Can NANOSEAL be applied during whitening treatment? 12

Other Q&A

- Q1. How many drops can be collected in one bottle? How many teeth can be treated with a single collection?...13
- Q2. What should be done for the maintenance of the attached dappen dish after use? 13
- Q3. Is there any dedicated tool available for application? 13
- Q4. How can the liquid stuck to a table or other dental equipment be removed?..... 13
- Q5. How can the liquid adhered to clothes be removed?..... 14
- Q6. How should NANOSEAL be stored?..... 14
- Q7. Can the once-frozen product still be used for application? 14
- Q8. What should I do if the tip of the nozzle is clogged?..... 14

Q&A about Product Characteristics

Q1. What type of dentin desensitizer is NANOSEAL?


A The products with the following three characteristics are primarily available as a dentin desensitizer¹⁾:

1 Nerve desensitization ... "Make the nerve less sensitive"

2 Coagulation ... "Coagulate protein to prevent movement of tissue fluids in the dentinal tubules."

3 Sealing ... "Seal the dentinal tubule."

Among them, **NANOSEAL is a product that seals the dentinal tubule.**




[Figure by the courtesy of Prof. Morioki Fujitani and Prof. Akira Senda (Operative Dentistry Course, School of Dentistry, Aichi Gakuin University)]

Q2. What are the new characteristics of NANOSEAL?

- A**
- ① **Immediately reduce hypersensitivity just after the application.**
 - ② **Protect from demineralization**
 - ③ **Promote remineralization**
 - ④ **Easier than ever! There is no need for drying, rubbing or light-curing before/after the application.**

Materials with
all of these characteristics
**have never been used
in traditional
dentin desensitizers!**



Q3. What is the composition and pH of NANOSEAL? Does it smell?


- A** NANOSEAL consists of two liquids, Liquid A and Liquid B.
Each composition is as follows:

	Composition	pH	odor
Liquid A	fluoroaluminosilicate glass dispersion liquid (paraben is used as a preservative), purified water	approx. 8.0-9.0	odorless
Liquid B	aqueous phosphate solution(approximately 10%), purified water	approx. 1.0	odorless
Mixture		approx. 2.0-3.5	odorless


Q&A about Product Characteristics

Q4. How is the acid-resistant nanoparticle layer formed?


A A solid nanoparticle layer incorporated with the tooth substance (thickness: approximately $1\ \mu\text{m}^2$) is formed after an immediate reaction with the tooth substance by the application of the mixture of Liquid A and Liquid B on the tooth surface.



《The nanoparticle layer seals the dental tubules³⁾》




Opened dental tubules



Sealed dental tubules after the application of NANOSEAL

Q5. What is the composition of the acid-resistant nanoparticle layer?

A It contains calcium fluoride, calcium phosphate, calcium silicate, and fluoroaluminosilicate glass, etc.




The nanoparticle layer contains fluoroaluminosilicate glass.

Q6. How can resistance to acids of the nanoparticle layer be provided?

A The acid resistance can be gained with the acid buffer capacity of fluoroaluminosilicate glass and calcium silicate.

Q7. How strong is the acid resistance of the nanoparticle layer?

A It has been reported that no change was observed in the acid resistance after 110-hour immersion in carbonated beverages, that is, cola drinks (pH 2.7).⁴⁾



It has been reported that NANOSEAL-treated tooth enamel and dentin exhibited superior acid resistance.

Q8. Does NANOSEAL contain fluorine?


A The fluoroaluminosilicate glass (nanoparticles) in Liquid A contains **fluorine almost equivalent to that present in typical acidulated phosphate fluoride (APF)** preparations, which implies that the nanoparticle layer also contains fluorine.

Q&A about Usage

Before application

Q1. How long can the mixture of Liquid A and Liquid B be used?

A Use it **within 1 minute** after mixing of the two liquids. As shown in the figure, the pH of the mixture increases over time, which reduces the reactivity. Therefore, the rate of sealing of the dentinal tubules may start to gradually decline 1 minute after the mixing, undermining the preventive effects against hypersensitivity.



Q2. Is there any problem if the collection ratio of Liquid A and Liquid B is not equal?

A If the difference is within the range of variation observed during one-drop collection, the volumes of the two liquids do not need to be exactly the same. However, a significant difference in the collection ratio (two drops or more vs. one drop) may attenuate the preventive effects against hypersensitivity.

Q3. Is any pre-application treatment required for the tooth surface other than cleaning?

A No treatment such as air-drying is required for the tooth surface. If cleaning the tooth surface is difficult and the surface is stained only with saliva or if there is a small amount of blood due to bleeding during scaling and root planning (SRP), NANOSEAL is still effective and the dentinal tubules can be sealed with the application without cleaning the surface.⁵⁾ However, we recommend removing massively accumulated plaque and tartar as much as possible.

During application

Q4. How should NANOSEAL be applied?

A Apply NANOSEAL sufficiently to the affected area. If hypersensitivity remains even after applying and rinsing the product with water, the area might have been left partially uncoated. In this case, apply the product again using a new mixture.



Q5. Should the surface be recoated?


- A** There is **no need** for recoating.
We rather recommend applying NANOSEAL sufficiently to the affected area to prevent uneven coating or partial coating.

Q6. Is there any problem if NANOSEAL contacts to the gingiva during application?

- A** **No problem, if NANOSEAL contacts to the gingiva.**
The pH of the mixture is approximately 2.0–3.5, which is equivalent to that of cola drinks, orange juice, and isotonic drinks.⁴⁾ (See the graph in Q1 on p. 7)

Q7. Do patients feel a bitter taste during application?

- A** **No, they don't feel a bitter taste.**
It may taste slightly sour, though.



Q8. Does the tooth hurt during application?

- A** **No, it does not.**
No irritation occurs because the acid-resistant nanoparticle layer is formed immediately after the application.

After application

Q9. Can the surface be rinsed with water immediately after the application?

- A** **The tooth can be rinsed with water soon after the application** because NANOSEAL immediately works on the tooth substance. Although it may taste slightly sour if left for a while after the application, it does not have adverse effects on the tooth substance and gingiva.


Q&A about Desensitizing Effect

Q1. What desensitizing effects of NANOSEAL have been reported against hypersensitivity?

A It has been reported that a **single application of NANOSEAL resulted in a significant response** (felt no pain at all) **in approximately 77% of the patients** and that these desensitizing effects were maintained even 2 weeks after the application, without showing a significant reduction, as depicted in the figure below.⁴⁾ Furthermore, since the nanoparticle layer is incorporated with the tooth substance, the applied NANOSEAL **will not easily come off**. Moreover, **the effects will not be attenuated**.

《Efficacy and durability of the desensitizing effects (in 30 patients)》

Terata R.:The Nippon Dental Review, 72(3), 89-94, 2012.



[Evaluation Procedure] Evaluation was performed based on pain that occurred due to 3-minute cold air stimulation generated by an air syringe placed approximately 1 cm away from the affected area or tactile stimulation using the tip of the dental probe. The efficacy was assessed using pain scores.

Q2. Will additional application be beneficial if no benefit has been observed with a single application?

A Improvements in the symptoms by additional application were reported in patients whose symptoms had been persistent at 1 or 2 weeks after the previous application of NANOSEAL.⁶⁾

《Change in the significant response rate》


Just after treatment	Two weeks after treatment
68.4% (in 65 of 95 patients)	96.8% (in 92 of 95 patients)

additionally applied

Q3. I am afraid that the nanoparticle layer may be removed by brushing.

A It cannot be removed so easily.

The presence of the nanoparticle layer was confirmed on the tooth surface even after 5-minute brushing (toothbrush pressure, 250 g; scrubbing method) using a toothpaste containing an abrasive.⁷⁾ However, excessively strong brushing may shave and wear away the tooth. Guidance for brushing should be provided for a patient, if necessary.



Q4. Can NANOSEAL be applied to the root surface with the remaining cementum?

A Yes, it can be applied.

Desensitizing effect can still be expected for the root surface with the remaining cementum as NANOSEAL works on the cementum and forms the nanoparticle layer.

Q5. Does NANOSEAL have the same effects if used on the surface previously coated with a resin-based dentin desensitizer?

A Effects cannot be expected because NANOSEAL does not act on resin.

We recommend removing the resin to apply the product to the dentinal surface.

Q&A about Application to the Prepared to

Q1. Can NANOSEAL be used for the prepared tooth?

A Yes, it can.

Apply NANOSEAL sufficiently to prevent uneven coating or partial coating.



Q2. Will the bond strength of resin cement and bonding materials be decreased on the prepared tooth surface coated with NANOSEAL?

A No, it will not.

It has been reported that application of NANOSEAL does not interfere with the bond strength of various resin cement and bonding materials used for the dentin.⁹⁾

Q3. Is phosphoric acid etching on the NANOSEAL-coated surface (enamel) effective?

A Yes, it is effective.

It has been demonstrated that phosphoric acid etching on the NANOSEAL-coated surface (enamel) increases the bond strength of the adhesive resin.

【Test Method】

A flat surface of the bovine tooth enamel where a masking tape is attached (thickness, 50 μ m; pore size, 3 mm) was used for the test. The inner surface was treated with phosphoric acid etching. The bond strength at a crosshead speed of 1 mm/min was evaluated using treated and untreated materials where Clearfil Mega Bond (Kuraray Noritake Dental Inc. (former Kuraray Medical Inc.)) was applied.

【Test Results】


Without phosphoric acid etching: approximately 10 Mpa/ With phosphoric acid etching: approximately 25 Mpa (Laboratory data provided by Nippon Shika Yakuhin Co., Ltd.)

Q4. Will the fitting of dental prostheses be affected by NANOSEAL?

A No, it will not.

The thickness of the nanoparticle layer is approximately 1 μm , which is much smaller compared to that required for dental prostheses (several tens of μm). The thickness will increase to only about 2 μm if the application and rinsing processes are repeated four times.¹⁰⁾

Repeated application
will not significantly
increase the thickness.

**Q5. Can temporary seal and cementation materials be used after the application of NANOSEAL?**


A All types of materials can be used.

No separating agent is required to use them.

Q6. Can NANOSEAL be applied during whitening treatment?

A Yes, it can.

It has been confirmed that NANOSEAL does not affect the whitening results¹¹⁾ and can be used both before and during whitening treatment. Furthermore, no coloring will occur after the application.



Other Q&A

Q1. How many drops can be collected in one bottle? How many teeth can be treated with a single collection?


A Approximately **100 drops** can be collected from Liquid A and Liquid B, respectively.

At least **two to three teeth** can be treated with one collection (one drop each).

Q2. What should be done for the maintenance of the attached dappen dish after use?

A Follow the procedure below;

◆ Dappen dish



Cleaning

After rinsing with water, wipe the dappen dish with a clean alcohol cotton ball.

Sterilization

Chemical sterilization with glutaral can be applied; however, autoclave sterilization should be avoided, as it may result in deformation of the dappen dish due to heat.

Q3. Is there any dedicated tool available for application?

A No dedicated tools are available. Use a cotton ball or brush for application. Use a cotton ball for easy application of a large amount to multiple teeth.

Q4. How can the liquid stuck to a table or other dental equipment be removed?

A Wipe off the liquid using a dampen cloth or a paper towel several times. Fluoroaluminosilicate glass used for this product consists of fine particles; therefore, white spots may still remain after wiping off the liquid. They can be removed by repeated wiping.


Q5. How can the liquid adhered to clothes be removed?

- A** Wipe off the liquid adhered to the clothes by getting them wet with water using a dampened cloth or a wet paper towel.
If possible, wash the adhered part by soaking it in water.

Q6. How should NANOSEAL be stored?

- A** Avoid direct sunlight and freezing. Store NANOSEAL at room temperature (1°C–30°C). Storage in a refrigerator is recommended if the temperature is expected to be 0°C or below.


NANOSEAL
gets frozen
at 0°C or below.



Q7. Can the once-frozen product still be used for application?

- A** Yes, it can.
However, coagulated glass components of Liquid A are slightly likely to form precipitates after being frozen. Shake it well before use to ensure even mixing.

Shake
Liquid A
bottle
well.




Q8. What should I do if the tip of the nozzle is clogged?

- A** Tighten the cap and place the bottle in the basket of an ultrasonic cleaner to run the cleaner for about 20 seconds. Ultrasonic vibration will clear the clog.

nanoseal User Guide

1 Dispense & Mix




Shake the Liquid A bottle gently.
Place an equal amount of Liquid A and B
in the dappen dish to mix them.

Apply the mixture
within one(1) minute after mixing.

★ In Liquid A dispensed after shaking the bottle gently, aggregated glass particles may rarely be found, which would not affect the properties of this product.


2 Apply



After cleaning tooth surface,
apply sufficiently.
(Repeat application for five (5)
to twenty(20) secs.)

Follow the same procedures
when applied to prepared tooth surfaces.

3 Rinse



Wash or rinse.

★ Please read "Instructions for use" carefully before use.

Literature

- 1) Fujitani Morioki, Senda Akira: Hypersensitive Dentin 2nd edition, Perfect Treatment Guide, Ishiyaku Pub, Inc., Tokyo, 2013. (in Japanese)
- 2) HAN Linlin, ISHIZAKI Hiroko, FUKUSHIMA Masayoshi, and OKIJI Takashi: Effect of a Prototype Fluoride-containing Tooth Surface Coating Material on Enamel and Dentin Surfaces; The Japanese Journal of Conservative Dentistry Vol.55, No.1 53-59, 2012
- 3) KASHIWAMURA Haruko, YANAGIDA Kenichi, NOMURA Yuko, OZAKI Masao, MOTOKAWA Wataru: Study on Dentin Tubules Occlusion of "a prototype fluoride-containing nanoseal material" ; The Japanese Journal of Pediatric Dentistry, Vol. 48, No.2:308, 2010. (in Japanese)
- 4) TERATA Rintaro: Nanoseal, A New Dental Material Forming an Acid-resistant Layer on Tooth Surfaces-the Desensitizing Effects and Potentials; The Nippon Dental Review, 72(3), 89-94, 2012 (in Japanese)
- 5) YAMAGUCHI Hiroyasu, YAHAGI Hozumi, ASANO Souei, YOKOTA Kazuyoshi and TSUNEKAWA Masayoshi: *In-vitro* Dentinal Tubule Occlusion and Clinical Desensitization by Application of "Nanoseal," The Japanese Journal of Conservative Dentistry, Vol.57, No.4 333-342, 2014
- 6) TERATA Rintaro: Clinical Evaluation of Dentin Desensitizer NANOSEAL; Program and Abstracts, the 140th Meeting of the Japanese Society of Conservative Dentistry, 191, 2014. (in Japanese)
- 7) NOMURA Yuko, KASHIWAMURA Haruko, YANAGIDA Kenichi, OZAKI Masao: Study on Caries Prevention of "a prototype fluoride-containing nanoseal material, -Acid Resistance and Abrasion Resistance when Dentin tubules are occluded; The Japanese Journal of Pediatric Dentistry, vol. 48 No.5: 595, 2010. (in Japanese)
- 8) NISHIKAWA Yoshiaki, KUWATA Masahiro, ed.: "Desensitizer NANOSEAL" ; Dental Diamond, 37(4), 137-141, 2012. (in Japanese)
- 9) OKADA Nobuo, SHIGA Hanae, CHIBA Fumiko, NODA Mamoru, YOKOTA Kazuyoshi, TSUNEKAWA Masayoshi.; Program and Abstracts, the 139th Meeting of the Japanese Society of Conservative Dentistry, 166, 2013. (in Japanese)
- 10) Han Rinrin, Fukushima Masayoshi, Okiji Takashi: Effect of a prototype fluoride-containing sealing material on the surface ultrastructure of enamel and dentin, The Japanese Journal of Conservative Dentistry, Vol.54, 126, 2011.
- 11) Lu Ling, Yasuo Kenzo, Onda Kohei, Yoshikawa Kazushi, Yamamoto Kazuyo : Influence of various dentin desensitizers on the effect of tooth whitening : Journal of Osaka Dental University 49(1), 49-60, 2015



【Packaging】 Liquid A 5mL, Liquid B 5mL, dappen dish 1 piece



Nippon Shika Yakuhin Co., Ltd.

4-7-24 Takezaki-cho, Shimonoseki-shi, Yamaguchi 750-0025 Japan

TEL : (81) 83-222-2221 FAX: (81) 83-222-2220

URL: <https://www.nishika.co.jp/english/>